

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human BMF / Bcl2 modifying factor Protein (His tag)

产品货号: TP06923

第三版

别名:BMF

描述:BMF(Bcl2 modifying factor) belongs to the BCL2 protein family. BCL2 family members form hetero- or homodimers and act as anti- or pro-apoptotic regulators that are involved in a wide variety of cellular activities. BMF contains a single BCL2 homology domain 3 (BH3), and has been shown to bind BCL2 proteins and function as an apoptotic activator. BMF is found to be sequestered to myosin V motors by its association with dynein light chain 2, which may be important for sensing intracellular damage and triggering apoptosis. Alternatively spliced transcript variants encoding different isoforms have been identified.

配方:PBS

分子量:22.3 kDa

序列:Met 1-Arg184

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.