



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human TLE3 Protein (aa 484-772, GST tag)

产品货号: TP06924

第三版

别名:ESG;ESG3;GRG3;HsT18976

描述:Immunoglobulin superfamily member 11(IGSF11) are expressed on the plasma membrane in the testis and brain. These IGSF proteins undergo final modifications during capacitation and/or the acrosome reaction. IGSF proteins share significant homology with endothelial cell-selective adhesion molecule and coxsackievirus and adenovirus receptor, which mediates cell attachment and homotypic intercellular interactions. In clinical, the IGSF11 has been reported to over expressed in colorectal cancers and hepatocellular carcinomas as well as intestinal-type gastric cancers compared to their corresponding non- cancerous tissues. The IGSF11 has also been found expressed abundantly in testis and ovary and the IGSF11 can be used as a candidate of cancer-testis antigen.

配方:PBS

分子量:58.5 kDa

序列:Ser 484-Tyr 772

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.