
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Rat CD133 / PROM1 / Prominin 1 Protein (Fc tag)

产品货号: TP06931

第三版

别名:PROM1

描述:CD133, also known as PROM1 and Prominin 1, is a pentaspan transmembrane glycoprotein which belongs to the prominin family. It localizes to membrane protrusions and is often expressed on adult stem cells. CD133 is known to play a role in maintaining stem cell properties by suppressing differentiation. CD133 binds cholesterol in cholesterol-containing plasma membrane microdomains. It is proposed to play a role in apical plasma membrane organization of epithelial cells. CD133 is also involved in regulation of MAPK and Akt signaling pathways. Mutations in PROM1 gene have been shown to result in retinitis pigmentosa and Stargardt disease. PROM1 gene is expressed from at least five alternative promoters that are expressed in a tissue-dependent manner. Expression of this gene is also associated with several types of cancer.

配方:PBS

分子量:55.1 kDa

序列:Asn171-Tyr424

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.