

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Recombinant Human BLNK / Ly-57 / SLP-65 Protein (His tag)

产品货号: TP06959

第三版

别名:AGM4;BASH;bca;BLNK-S;LY57;SLP-65;SLP65

**描述:**B-cell linker protein, also known as B-cell adapter containing a SH2 domain protein, B-cell adapter containing a Src homology 2 domain protein, Cytoplasmic adapter protein, Src homology 2 domain-containing leukocyte protein of 65 kDa, SLP-65 and BLNK, is a cytoplasm and cell membrane protein which contains one SH2 domain. BLNK is expressed in B-cell lineage and fibroblast cell lines. Highest levels of expression is in the spleen, with lower levels in the liver, kidney, pancreas, small intestines and colon. BLNK functions as a central linker protein that bridges kinases associated with the B-cell receptor (BCR) with a multitude of signaling pathways, regulating biological outcomes of B-cell function and development. BLNK plays a role in the activation of ERK / EPHB2, MAP kinase p38 and JNK. BLNK modulates AP1 activation. It is important for the activation of NF-kappa-B and NFAT. BLNK plays an important role in BCR-mediated PLCG1 and PLCG2 activation and Ca<sup>2+</sup> mobilization and is required for trafficking of the BCR to late endosomes. BLNK may be required for the RAC1-JNK pathway. It plays a critical role in orchestrating the pro-B cell to pre-B cell transition. BLNK also plays an important role in BCR-induced B-cell apoptosis. Defects in BLNK are the cause of agammaglobulinemia type 4 (AGM4) which is a primary immunodeficiency characterized by profoundly low or absent serum antibodies and low or absent circulating B cells due to an early block of B-cell development.

**配方:**PBS

**分子量:**53 kDa

**序列:**Met 1-Ser 456

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.