

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Recombinant Human OMGP / OMG Protein (aa 1-416, His tag)

产品货号: TP06960

第三版

别名:OMGP

**描述:**Oligodendrocyte-myelin glycoprotein, also known as OMG and OMGP, is a cell membrane protein which contains eight LRR (leucine-rich) repeats. OMG / OMGP is a glycosylphosphatidylinositol-anchored protein expressed by neurons and oligodendrocytes in the central nervous system (CNS). OMG / OMGP is a cell adhesion molecule contributing to the interactive process required for myelination in the central nervous system. OMG / OMGP play roles in both the developing and adult central nervous system. OMG / OMGP participats in growth cone collapse and inhibition of neurite outgrowth through its interaction with NgR, the receptor for Nogo. This function requires its leucine-rich repeat domain, a highly conserved region in OMgp during mammal evolution. OMG / OMGP leucine-rich repeat domain is also implicated in the inhibition of cell proliferation. OMG / OMGP may also be involved in the formation and maintenance of myelin sheaths. Cell proliferation, neuronal sprouting and myelination are crucial processes involved in brain development and regeneration after injury.

**配方:**PBS

**分子量:**46 kDa

**序列:**Met 1-Pro 416

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.