



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human SerpinB1 / ELANH2 Protein (His tag)

产品货号: TP06961

第三版

别名:EI;ELANH2;HEL-S-27;HEL57;LEI;M/NEI;MNEI;PI-2;PI2

**描述:**SerpinB1 is an endogenous inhibitor of serine proteases recognized for its anti-inflammatory and host-protective properties. Serum SerpinB1 levels are elevated in patients with type 2 diabetes compared with that in healthy subjects and are negatively correlated with serum LDL-C. serpinB1 as an endogenous protein that can potentially be harnessed to enhance functional beta cell mass in patients with diabetes.

**配方:**PBS

**分子量:**44 kDa

**序列:**Met 1-Pro 379

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.