

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse Smad5 Protein (His & GST tag)

产品货号: TP06982

第三版

别名:1110051M15Rik;AI451355;Dwf-C;Madh5;MusMLP

描述:SMAD5 is a member of the SMAD family. Members of this family mediate signal transduction by the TGF-beta/activin/BMP-2/4 cytokine superfamily from receptor Ser/Thr protein kinases at the cell surface to the nucleus. SMAD5 is involved in the TGF-beta signaling pathway that results in an inhibition of the proliferation of hematopoietic progenitor cells. It is also involved in cell signalling and modulates signals of bone morphogenetic proteins (BMP's). The binding of ligands causes the oligomerization and phosphorylation of the SMAD5 protein. SMAD5 is a receptor regulated SMAD (R-SMAD) and is activated by bone morphogenetic protein type 1 receptor kinase.

配方:PBS

分子量:79.9 kDa

序列:Thr2-Ser465

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.