

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Calumenin Protein (His Tag)

产品货号: TP06988

第三版

别名:

描述: Calumenin belongs to the CREC family. It contains 6 EF-hand domains. Calumenin is expressed in skeletal muscle (at protein level). Calumenin interacts with GGCX and RYR1 in the presence of calcium ions, but not in the presence of EDTA. Calumenin is involved in regulation of vitamin K-dependent carboxylation of multiple N-terminal glutamate residues. It seems to inhibit gamma-carboxylase GGCX. Calumenin also binds 7 calcium ions with a low affinity and may modulate calcium release from the sarcoplasmic reticulum.

配方: PBS

分子量: 35.9 kDa

序列: Met 1-Arg311

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.