



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Rat CLEC4B2 / mDCAR1 Protein (His tag)

产品货号: TP06990

第三版

别名:CLEC4B2;Aplra1

描述:Clec4b2, also known as mDCAR1, is a member of the DCIR/DCAR family. Expression of Clec4b2 was strongly tissue dependent. Clec4b2 expression on DCs was restricted to the CD8(+) DC subset in spleen and thymus and on subpopulations of CD11b(+) myeloid cells in bone marrow and spleen, whereas the molecule was not detectable on both cell types in lymph nodes and peripheral blood. Clec4b2 is a functional receptor on cells of the immune system and provides further insights into the regulation of immune responses by CLR.

配方:PBS

分子量:21.8 kDa

序列:Met44-Leu208

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.