

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human FSTL5 Protein (His tag)

产品货号: TP07009

第三版

别名:FSTL5;KIAA1263

描述:FSTL5 may have molecular function (calcium ion binding) and to localize in various compartments (cytoplasm, extracellular space, extracellular region). FSTL5 expression denoted a dismal prognosis both within and across medulloblastoma subgroups. FSTL5 gene is well expressed, 1.0 times the average gene in this release. The sequence of this gene is defined by 120 GenBank accessions from 113 cDNA clones, some from brain, cerebellum, eye, melanotic melanoma, skin, amygdala, breast and 24 other tissues. FSTL5 gene contains 27 distinct introns. The addition of FSTL5 immunohistochemistry to existing molecular stratification schemes constitutes a reliable and cost-effective tool for prognostication in future clinical trials of medulloblastoma.

配方:PBS

分子量:94.9 kDa

序列:Gln32-Ala847

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.