



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human LRG1 Protein (Fc tag)

产品货号: TP07022

第三版

别名:HMFT1766;LRG

描述:LRG1 belongs to the leucine-rich repeat (LRR) family. Members of this family are involved in protein-protein interaction, signal transduction, and cell adhesion and development. LRG1 is expressed during granulocyte differentiation. It contains 4 LIM zinc-binding domains and 1 Rho-GAP domain.

配方:PBS

分子量:61.3 kDa

序列:Met 1-Gln347

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.