
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human SNAP25 / SUP Protein (His tag)

产品货号: TP07057

第三版

别名:bA416N4.2;dJ1068F16.2;RIC-4;RIC4;SEC9;SNAP;SNAP-25

描述: Synaptosomal-associated protein 25, also known as Super protein, Synaptosomal-associated 25 kDa protein, SNAP25 and SNAP, is a cytoplasm and cell membrane protein which belongs to the SNAP-25 family. SNAP25 / SUP contains 2 t-SNARE coiled-coil homology domains. SNAP25 / SUP is a membrane bound protein anchored to the cytosolic face of membranes via palmitoyl side chains in the middle of the molecule. SNAP25 / SUP protein is a component of the SNARE complex, which is proposed to account for the specificity of membrane fusion and to directly execute fusion by forming a tight complex that brings the synaptic vesicle and plasma membranes together. SNAP25 / SUP is a Q-SNARE protein contributing two α -helices in the formation of the exocytotic fusion complex in neurons where it assembles with syntaxin-1 and synaptobrevin. SNAP25 / SUP is involved in the molecular regulation of neurotransmitter release. It may play an important role in the synaptic function of specific neuronal systems. SNAP25 / SUP associates with proteins involved in vesicle docking and membrane fusion. SNAP25 / SUP regulates plasma membrane recycling through its interaction with CENPF. SNAP25 / SUP inhibits P/Q- and L-type voltage-gated calcium channels located presynaptically and interacts with the synaptotagmin C2B domain in Ca^{2+} -independent fashion. In glutamatergic synapses SNAP25 / SUP decreases the Ca^{2+} responsiveness, while it is naturally absent in GABAergic synapses.

配方: PBS

分子量: 24.8 kDa

序列: Met 1-Gly 206

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.