
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse Peroxiredoxin 1 / PRDX1 Protein (His tag)

产品货号: TP07061

第三版

别名: MSP23; NkefA; OSF-3; OSF3; PAG; Paga; Prdx1; prx1; Prx1; Tdpx2; TDX2; TPxA

描述: Peroxiredoxin-1, also known as Thioredoxin peroxidase 2, Natural killer cell-enhancing factor A, PRDX1, and PAGA, is a member of the ahpC/TSA family. Peroxiredoxin-1 is constitutively expressed in most human cells. It is induced to higher levels upon serum stimulation in untransformed and transformed cells. Peroxiredoxins (PRDXs) are a family of antioxidant enzymes that are also known as scavengers of peroxide in mammalian cells. The overexpression of Peroxiredoxin-1, which is one of the peroxiredoxins that is a ubiquitously expressed protein, was related to a poor prognosis in several types of human cancers. Peroxiredoxin-1 is involved in redox regulation of the cell. It reduces peroxides with reducing equivalents provided through the thioredoxin system but not from glutaredoxin and may play an important role in eliminating peroxides generated during metabolism. Peroxiredoxin-1 might participate in the signaling cascades of growth factors and tumor necrosis factor- α by regulating the intracellular concentrations of H₂O₂. The reduced Peroxiredoxin-1 expression is an important factor in esophageal squamous cancer progression and could serve as a useful prognostic marker.

配方: PBS

分子量: 23.5 kDa

序列: Met 1-Lys 199

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 μ g of protein (determined by LAL method)

存储: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.