

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human CFHR1 Protein (His Tag)

产品货号: TP07074

第三版

别名:CFHL;CFHL1;CFHL1P;CFHR1;CFHR1P;FHR1;H36-1;H36-2;HFL1;HFL2

描述:CFHR1 is a secreted protein belonging to the complement factor H protein family. The human complement factor H protein family consists of the complement and immune regulators factor H, the factor H-like protein 1 (FHL-1) and five factor H-related proteins (CFHR-1 to -5). Members of the H-related protein family are exclusively composed of individually folded protein domains, termed short consensus repeats (SCRs) or complement control modules. CFHR1 binds to *Pseudomonas aeruginosa* elongation factor Tuf together with plasminogen, which is proteolytically activated. CFHR1 might be involved in complement regulation. It can associate with lipoproteins and may play a role in lipid metabolism.

配方:PBS

分子量:37.2 kDa

序列:Met 1-Arg330

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.