
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human S100A7 / PSOR1 Protein

产品货号: TP07106

第三版

别名:PSOR1;S100A7c

描述:Protein S100-A7, also known as S100 calcium-binding protein A7, Psoriasin, S100A7, and PSOR1, is a secreted protein which belongs to the S-100 family. S100A7 was first isolated from skin involved by psoriasis, which can be induced in cultured squamous epithelial cells. S100A7 is expressed by both normal cultured and malignant keratinocytes and malignant breast epithelial cells within ductal carcinoma in situ, suggesting an association with abnormal pathways of differentiation. S100A7 plays a role in the pathogenesis of inflammatory skin disease, as a chemotactic factor for hematopoietic cells. It also plays a role in early stages of breast tumor progression in association with the development of the invasive phenotype.

配方:PBS

分子量:11.6 kDa

序列:Met 1-Gln 101

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.