
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human AMIGO2 Protein (His tag)

产品货号: TP07121

第三版

别名:ALI1;AMIGO-2;DEGA

描述:AMIGO2 contains Ig-like C2-type (immunoglobulin-like) domain, 6 LRR (leucine-rich) repeats, 1 LRRCT domain and 1 LRRNT domain. It belongs to the immunoglobulin superfamily, AMIGO family. AMIGO2 may mediate homophilic as well as heterophilic cell-cell interaction with AMIGO1 or AMIGO3. It is required for depolarization-dependent survival of cultured cerebellar granule neurons. AMIGO2 may also contribute to signal transduction through its intracellular domain. It may play a role in the tumorigenesis of a subset of gastric adenocarcinomas. AMIGO2 is highly expressed in breast, ovary, cervix, and uterus.

配方:PBS

分子量:42.3 kDa

序列:Met 1-Thr398

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.