
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human ECE1 Protein (His tag)

产品货号: TP07180

第三版

别名: ECE; RP3-329E20.1

描述: Endothelin-converting enzyme 1, also known as ECE-1, is a single-pass type I I membrane protein which belongs to the peptidase M13 family. ECE-1 converts big endothelin-1 to endothelin-1. ECE-1 is a membrane metalloprotease that generates endothelin from its direct precursor big endothelin. Four isoforms of ECE-1 are produced from a single gene through the use of alternate promoters. These isoforms share the same extracellular catalytic domain and contain unique cytosolic tails, which results in their specific subcellular targeting. All isoforms of ECE-1 are expressed in umbilical vein endothelial cells, polynuclear neutrophils, fibroblasts, atrium cardiomyocytes and ventricles. Isoforms A, B and C of ECE-1 are also expressed in placenta, lung, heart, adrenal gland and pheochromocytoma; isoforms A and C of ECE-1 in liver, testis and small intestine; isoform B, C and D of ECE-1 in endothelial cells and umbilical vein smooth muscle cells; isoforms C and D in saphenous vein cells, and isoform C in kidney. Defects in ECE1 are a cause of Hirschsprung disease, cardiac defects and autonomic dysfunction. It is a form of Hirschsprung disease with skip-lesions defects, craniofacial abnormalities and other dysmorphic features, and autonomic dysfunction.

配方: PBS**分子量:** 80 kDa**序列:** Gln 90-Trp 770**纯度:** > 95% by HPLC**浓度:****内毒素:** <1.0 EU per 1 ug of protein (determined by LAL method)**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.