
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Rat CLEC4F / CLECSF13 Protein

产品货号: TP07247

第三版

别名:CLEC4F;Clecsf13;Kclr

描述:CLEC4F, a member of C-type lectins, was firstly purified from rat liver extract with high binding affinity to fucose, galactose and N-acetylgalactosamine, and un-sialylated glucosphingolipids with GalNAc or Gal terminus. However, the biological functions of CLEC4F have not been elucidated. Histochemical staining showed that mouse CLEC4F $\dot{\bar{m}}$ $\dot{\bar{m}}$ CLEC4F $\dot{\bar{m}}$ $\dot{\bar{m}}$ is only expressed on F4/80+ cells localized in liver, and is undetectable in bone marrow, spleen, lymph nodes, or other tissues in adult mice. However, mCLEC4F is detected in the liver of embryonic day 11.5 (E11.5), which is 1.5 day earlier than the formation of liver (E10) and is 3.5 day earlier than the formation of bone marrow (E15-16). Moreover, recombinant mCLEC4F.Fc binds to alpha-galactoceramide in a Ca⁺⁺-dependent manner, and both galactose and ceramide can partially inhibit CLEC4F.Fc binding to alpha-galactoceramide. Interestingly, mCLEC4F-deficient (mCEC4F k/o) mice produced far less cytokines than wild type littermates after intravenous injection of alpha-galactoceramide. This suggests that mCLEC4F is not only a specific marker for Kupffer cells, but is also critical for the presentation of glycolipid antigen to NKT cells.

配方:PBS**分子量:**53.7 kDa**序列:**Arg70-Ser550**纯度:**> 95% by HPLC**浓度:****内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.