
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human 14-3-3 tau / 14-3-3 theta / YWHAQ Protein (GST tag)

产品货号: TP07252

第三版

别名:1C5;41701;HS1

描述:G protein-coupled receptor kinase 5, also known as G protein-coupled receptor kinase GRK5 and GRK5, is a member of the protein kinase superfamily, AGC Ser/Thr protein kinase family and GPRK subfamily. GRKs specifically phosphorylate agonist-occupied G protein-coupled receptors at the inner surface of the plasma membrane (PM), leading to receptor desensitization. GRKs utilize a variety of mechanisms to bind tightly, and sometimes reversibly, to cellular membranes. GRKs play an important role in mediating agonist-specific desensitization of numerous G protein-coupled receptors. GRK5 contains one AGC-kinase C-terminal domain, one protein kinase domain and one RGS domain. GRK5 specifically phosphorylates the activated forms of G protein-coupled receptors. Phospholipid-stimulated autophosphorylation may represent a novel mechanism for membrane association and regulation of GRK5 activity. GRK5 deficiency significantly exaggerates microgliosis and astrogliosis in the presence of an inflammatory initiator, such as the excess fibrillar Abeta and the subsequent active inflammatory reactions. GRK5 deficiency has been linked to early Alzheimer's disease in humans and mouse models of the disease.

配方:PBS

分子量:54.6 kDa

序列:Met 1-Asn 245

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.