

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human ALK-1 / ACVRL1 Protein (Fc tag)

产品货号: TP07267

第三版

别名:ACVRLK1;ALK-1;ALK1;HHT;HHT2;ORW2;SKR3;TSR-1

**描述:**Activin A receptor, type II-like 1 (ACVRL1), also known as ALK-1 (activin receptor-like kinase 1), is an endothelial-specific type I receptor of the TGF-beta (transforming growth factor beta) receptor family of ligands. On ligand binding, a heteromeric receptor complex forms consisting of two type II and two type I transmembrane serine/threonine kinases. ACVRL1 protein is expressed in certain blood vessels of kidney, spleen, heart and intestine, serving as an important role during vascular development. Mutations in ACVRL1 gene are associated with hemorrhagic telangiectasia type 2, also known as Rendu-Osler-Weber syndrome 2 and vascular disease.

**配方:**PBS

**分子量:**37.4 kDa

**序列:**Met 1-Gln 118

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.