
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Mouse MBL-2 / MBL Protein (His Tag)

产品货号: TP07292

第三版

别名:L-MBP;MBL;MBL-C;MBP-C

描述:MBL (mannose-binding lectin) is primarily a liver-derived collagen-like serum protein, which binds sugar structures on microorganisms and on dying host cells and is one of the four known mediators that initiate activation of the complement system via the lectin pathway. MBL and the ficolins (Ficolin-1, Ficolin-2 and Ficolin-3) are soluble collagen-like proteins that are involved in innate immune defence. They bind sugar structures or acetylated compounds present on microorganisms and on dying host cells and they initiate activation of the lectin complement pathway in varying degrees. MBL2 encodes the mannose-binding lectin, which is a key player in the innate immune system and has recently been found to play a role in development of type 1 diabetes and gestational diabetes mellitus. Common variant alleles situated both in promoter and structural regions of the MBL2 gene influence the stability and the serum concentration of the protein. Several polymorphisms in the promoter and structural regions of MBL2 adversely affect the plasma concentration and oligomeric state of MBL. The possession of mutant alleles has been linked to disease outcome for a variety of bacterial and viral infections. Mutant MBL2 haplotypes have been linked to disease progression and response to therapy in HCV infection.

配方:PBS

分子量:26.3 kDa

序列:Glu19-Asp244

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.