



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse Cathepsin H Protein (His Tag)

产品货号: TP07342

第三版

别名:AL022844 Protein, Mouse

描述:Proinflammatory cytokines decreased the expression of CTSH in human islets and primary rat beta-cells, and overexpression of CTSH protected insulin-secreting cells against cytokine-induced apoptosis.

配方:PBS

分子量:36.6 kDa

序列:Ala21-Val333

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.