
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human IL17RC Protein (aa 1-454, His tag)

产品货号: TP07347

第三版

别名:IL17-RL;IL17RL;UNQ6118/PRO20040/PRO38901

描述:The hypomethylation within the IL17RC gene promoter in peripheral blood is not suitable for use as a clinical biomarker of AMD. This study highlights the need for considerable replication of epigenetic association studies prior to clinical application. methylation of IL17RC could play as a marker in choroidal neovascularization (CNV) and degeneration of retinal pigment epithelium (RPE) cells in vitro.

配方:PBS

分子量:49.6 kDa

序列:Met 1-Ala 454

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.