

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human VEGFR3 / FLT4 Protein (Fc tag)

产品货号: TP07349

第三版

别名:FLT-4;FLT-41;FLT41;LMPH1A;PCL;VEGF Receptor 3;VEGFR-3;VEGFR3

描述: Vascular endothelial growth factor receptor 3 (VEGFR3), also known as FLT-4, together with the other two members VEGFR1 (FLT-1) and VEGFR2 (KDR/Flk-1) are receptors for vascular endothelial growth factors (VEGF) and belong to the class III subfamily of receptor tyrosine kinases (RTKs). The VEGFR3 protein is expressed mainly on lymphatic vessels but it is also up-regulated in tumor angiogenesis. Mutations in VEGFR3 have been identified in patients with primary lymphoedema. The VEGF-C/VEGF-D/VEGFR3 signaling pathway may provide a target for antilymphangiogenic therapy in prostate cancer, breast cancer, gastric cancer, lung cancer, non-small cell lung cancer (NSCLC), and so on. Immune Checkpoint Immunotherapy Cancer Immunotherapy Targeted Therapy

配方: PBS

分子量: 111 kDa

序列: Met 1-Ile 776

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.