

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human KLK3 / PSA / Kallikrein-3 Protein (His tag)

产品货号: TP07379

第三版

别名:APS;HK3;KLK2A1;PSA

**描述:**KLK3(Kallikrein 3), also known as Prostate Specific Antigen (PSA), kallikrein-related peptidase 3, Gamma-seminoprotein, is a secreted protein of the glandular kallikrein subfamily of serine proteases.KLK3 contains one peptidase S1 domain.KLK3 is a glycoprotein produced almost exclusively by the prostate gland. Growing evidence suggests that many kallikreins are implicated in carcinogenesis and some have potential as novel cancer and other disease biomarkers.Immune Checkpoint Immunotherapy Cancer Immunotherapy Targeted Therapy

**配方:**PBS

**分子量:**28.3 kDa

**序列:**Met 1-Pro 261

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.