
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Cystatin SN / CST1 Protein (His tag)

产品货号: TP07382

第三版

别名:CST1

描述: Cystatin-SN, also known as Cystatin-SA-I, Cystatin-1, Salivary cystatin-SA-1 and CST1, is a secreted protein which belongs to the cystatin family. The cystatin superfamily encompasses proteins that contain multiple cystatin-like sequences. Some of the members are active cysteine protease inhibitors, while others have lost or perhaps never acquired this inhibitory activity. Cystatin-SN / CST1 is expressed in submandibular and sublingual saliva but not in parotid saliva (at protein level). Cystatin-SN / CST1 is also expressed in saliva, tears, urine and seminal fluid. Human saliva appears to contain several cysteine proteinase inhibitors that are immunologically related to cystatin S but that differ in their specificity due to amino acid sequence differences. Human salivary cystatins include Cystatin-S, Cystatin-S1, Cystatin-S2, Cystatin-SA, Cystatin-SN, Cystatin-C and Cystatin-D. Cystatin-SN is a much better inhibitor of papain and dipeptidyl peptidase I than is cystatin-S, although both inhibit ficin equally well.

配方: PBS

分子量: 15.6 kDa

序列: Met 1-Ser 141

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.