
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Cystatin SA / CST2 Protein (His tag)

产品货号: TP07388

第三版

别名:CST2;MGC71924

描述: Cystatin-SA, also known as Cystatin-2, Cystatin-S5 and CST2, is a secreted protein which belongs to the cystatin family. Cystatin-2 / CST2 is expressed in submandibular and sublingual saliva but not in parotid saliva (at protein level). It is also expressed in submandibular gland and parotid gland. The cystatin superfamily encompasses proteins that contain multiple cystatin-like sequences. Some of the members are active cysteine protease inhibitors, while others have lost or perhaps never acquired this inhibitory activity. There are three inhibitory families in the superfamily, including the type 1 cystatins (stefins), type 2 cystatins and the kininogens. The type 2 cystatin proteins are a class of cysteine proteinase inhibitors found in a variety of human fluids and secretions, where they appear to provide protective functions. The CST1, CST2, CST4, and CST5 are expressed in differential, tissue-specific patterns. Expression of CST2 and CST5 is restricted to the submandibular and parotid glands, while CST1 and CST4 are expressed in these tissues and in the lacrimal gland.

配方: PBS

分子量: 16 kDa

序列: Met 1-Ala 141

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.