

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Recombinant Mouse 14-3-3 sigma / Stratifin / YWHAS Protein

产品货号: TP07403

第三版

别名:Er;Mme1;Ywhas

**描述:**14-3-3 protein sigma (YWHAS), also known as stratifin (SFN) and epithelial cell marker protein 1, is a member of the 14-3-3 proteins which are a family of conserved regulatory molecules expressed in all eukaryotic cells. The name 14-3-3 refers to the particular elution and migration pattern of these proteins on DEAE-cellulose chromatography and starch-gel electrophoresis. The 14-3-3 proteins eluted in the 14th fraction of bovine brain homogenate and were found on positions 3.3 of subsequent electrophoresis. There are seven genes that encode 14-3-3s in most mammals. 14-3-3 proteins have been identified as adapter proteins implicated in the regulation of a large spectrum of both general and specialized signaling pathways. More than 100 signaling proteins have been reported as 14-3-3 ligands including kinases, phosphatases, and transmembrane receptors, and the binding generally results in the modulation of the activity of the binding partner. YWHAS exists as a homodimer and is present mainly in tissues enriched in stratified squamous keratinising epithelium. YWHAS has been reported to interact with KRT17 and GAB2, and may regulate protein synthesis and epithelial cell growth by stimulating Akt/mTOR pathway upon binding to KRT17. Additionally, YWHAS (SFN) may also act as a p53-regulated inhibitor of G2/M progression.

**配方:**PBS

**分子量:**27.7 kDa

**序列:**Met 1-Ser 248

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.