
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human KIR2DL4 / CD158D Protein

产品货号: TP07490

第三版

别名:CD158D;G9P;KIR;KIR-103AS;KIR103;KIR103AS;KIR2DL4

描述:KIR2DL4, also known as CD158d, is a member of the killer cell Ig-like receptor (KIR) family. KIRs are transmembrane glycoproteins expressed by natural killer cells and subsets of T cells. The KIR genes are polymorphic and highly homologous. KIR2DL4 is expressed in all NK cells and some T cells. KIR2DL4 activates the cytotoxicity of NK cells, despite the presence of an immunoreceptor tyrosine-based inhibition motif (ITIM) in its cytoplasmic tail. The ITIM was not necessary for activation of lysis by KIR2DL4. The activation signal of KIR2DL4 was sensitive to inhibition by another ITIM-containing receptor. The activation-deficient mutant of KIR2DL4 inhibited the signal delivered by the activating receptor CD16.

配方:PBS

分子量:25.1 kDa

序列:Met 1-His242

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.