

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human RTN4 / NOGO-A Protein (GST tag)

产品货号: TP07502

第三版

别名:ASY;DGU;DKFZp781L1143;HIGM4;Nbla00271;Nbla10545;NI220/250;NOGO;NOGO-A;Nogo-B;Nogo-C;NOGOC;NSP;NSP-CL;RTN-X;RTN4-A;RTN4-B1;RTN4-B2;RTN4-C;UDG;UNG1;UNG15;UNG2

**描述:**Reticulon-4, also known as Focin, Neurite outgrowth inhibitor, Nogo protein, Neuroendocrine-specific protein, Neuroendocrine-specific protein C homolog, RTN-x, Reticulon-5 and RTN4, is a multi-pass membrane protein which contains one reticulon domain. Isoform 1 of RTN4 is specifically expressed in brain and testis and weakly in heart and skeletal muscle. Isoform 2 of RTN4 is widely expressed except for the liver. Isoform 3 of RTN4 is expressed in brain, skeletal muscle and adipocytes. Isoform 4 of RTN4 is testis-specific. Reticulon-4 / RTN4 is a developmental neurite growth regulatory factor with a role as a negative regulator of axon-axon adhesion and growth, and as a facilitator of neurite branching. Reticulon-4 / RTN4 regulates neurite fasciculation, branching and extension in the developing nervous system. Reticulon-4 / RTN4 is involved in down-regulation of growth, stabilization of wiring and restriction of plasticity in the adult CNS. It regulates the radial migration of cortical neurons via an RTN4R-LINGO1 containing receptor complex. Isoform 2 of RTN4 reduces the anti-apoptotic activity of Bcl-xl and Bcl-2. This is likely consecutive to their change in subcellular location, from the mitochondria to the endoplasmic reticulum, after binding and sequestration. Isoform 2 and isoform 3 of RTN4 inhibit BACE1 activity and amyloid precursor protein processing.

**配方:**PBS

**分子量:**46.2 kDa

**序列:**Met 1-Val 185

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.