
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human SULT1E1 / ST1E1 Protein (His tag)

产品货号: TP07507

第三版

别名:EST;EST-1;MGC34459;ST1E1;STE

描述:Estrogen sulfotransferase, also known as Sulfotransferase, estrogen-preferring, Sulfotransferase 1E1, SULT1E1 and ST1E1, is a cytoplasm enzyme which belongs to the sulfotransferase 1 family. Sulfotransferase enzymes catalyze the sulfate conjugation of many hormones, neurotransmitters, drugs, and xenobiotic compounds. These cytosolic enzymes are different in their tissue distributions and substrate specificities. SULT1E1 may control the level of the estrogen receptor by sulfurylating free estradiol. SULT1E1 maximally sulfates beta-estradiol and estrone at concentrations of 20 nM. SULT1E1 also sulfates dehydroepiandrosterone, pregnenolone, ethinylestradiol, equalenin, diethylstilbesterol and 1-naphthol, at significantly higher concentrations; however, cortisol, testosterone and dopamine are not sulfated. SULT1E1 is a key enzyme in estrogen homeostasis. It plays a central role in the prevention and development of human disease.

配方:PBS

分子量:36 kDa

序列:Asn 2-Ile 294

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.