
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human PCSK1 / NEC1 Protein (His tag)

产品货号: TP07521

第三版

别名: BMIQ12; NEC1; PC1; PC3; SPC3

描述: Neuroendocrine convertase 1, also known as Prohormone convertase 1, Proprotein convertase 1, PCSK1 and NEC1, is an enzyme which belongs to the peptidase S8 family and Furin subfamily. PCSK1 is an enzyme that performs the proteolytic cleavage of prohormones to their intermediate (or sometimes completely cleaved) forms. It is present only in neuroendocrine cells such as brain, pituitary and adrenal, and most often cleaves after a pair of basic residues within prohormones but can occasionally cleave after a single arginine. It binds to a protein known as proSAAS, which also represents its endogenous inhibitor. PCSK1 is involved in the processing of hormone and other protein precursors at sites comprised of pairs of basic amino acid residues. PCSK1 substrates include POMC, renin, enkephalin, dynorphin, somatostatin and insulin. Defects in PCSK1 are the cause of proprotein convertase 1 deficiency (PC1 deficiency). PC1 deficiency is characterized by obesity, hypogonadism, hypoadrenalism, reactive hypoglycemia as well as marked small-intestinal absorptive dysfunction. It is due to impaired processing of prohormones.

配方: PBS

分子量: 57.4 kDa

序列: Met 1-Arg 617

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.