
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human NMNAT2 / NMNAT-2 Protein (His tag)

产品货号: TP07523

第三版

别名:C1orf15;PNAT2

描述:NMNAT2, also known as NMNAT-2, belongs to the nicotinamide mononucleotide adenylyltransferase (NMNAT) enzyme family. NMNAT is a central enzyme in NAD⁺ biosynthesis, transferring the adenylyl moiety of ATP to nicotinamide mononucleotide (NMN) or nicotinic acid mononucleotide (NaMN) resulting in the formation of NAD⁺ or NaAD⁺ and the release of pyrophosphate. NMNAT2 is predominantly expressed in human pancreas, insulinoma as well as in the brain, especially in the cerebrum, cerebellum, occipital lobe, frontal lobe, temporal lobe and putamen. Immunofluorescence microscopy localized endogenous NMNAT2 to the Golgi apparatus in human cell line. Endogenous NMNAT2 seem to be a labile axon survival factor, because specific depletion of NMNAT2 is sufficient to induce Wallerian-like degeneration of uninjured axons which endogenous NMNAT1 and NMNAT3 cannot prevent. Thus endogenous NMNAT2 represents an exciting new therapeutic target for axonal disorders.

配方:PBS

分子量:35.8 kDa

序列:Met 1-Gly307

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.