

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human DUSP3 / VHR Protein

产品货号: TP07569

第三版

别名:VHR

**描述:** Vaccinia H1-related phosphatase (VHR) is classified as a dual-specificity phosphatase (DUSP), and the other name is dual-specificity phosphatase 3 (DUSP3). DUSPs are a heterogeneous group of protein phosphatases that can dephosphorylate both phosphotyrosine and phosphoserine/phosphothreonine residues within the one substrate. Unlike typical DUSPs, VHR lacks mitogen-activated protein kinase (MAPK)-binding domain, and shows poor activity against MAPKs. VHR often act on bisphosphorylated protein substrates, it displays a strong preference for dephosphorylating phosphotyrosine residues over phosphothreonine residues. VHR has been identified as a novel regulator of extracellular regulated kinases (ERKs). VHR is responsible for the rapid inactivation of ERK following stimulation and for its repression in quiescent cells. VHR is a negative regulator of the Erk and Jnk pathways in T cells and, therefore, may play a role in aspects of T lymphocyte physiology that depend on these kinases.

**配方:**PBS

**分子量:**20.6 kDa

**序列:**Met 1-Pro 185

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.