

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse S100A5 Protein (His tag)

产品货号: TP07638

第三版

别名:S100D9

描述:S100 protein is a family of low molecular weight protein found in vertebrates characterized by two EF-hand calcium-binding motifs. There are at least 21 different S100 proteins, and the name is derived from the fact that the protein is 100% soluble in ammonium sulfate at neutral pH. Most S100 proteins are disulfide-linked homodimer, and is normally present in cells derived from the neural crest, chondrocytes, macrophages, dendritic cells, etc. S100 proteins have been implicated in a variety of intracellular and extracellular functions. They are involved in regulation of protein phosphorylation, transcription factors, the dynamics of cytoskeleton constituents, enzyme activities, cell growth and differentiation, and the inflammatory response. Protein S100-A5, also known as Protein S-100D, S100 calcium-binding protein A5, S100A5 and S100D, is a member of the S100 family which contains two EF-hand domains. S100A5 is also a novel member of the EF-hand superfamily of calcium-binding proteins that is poorly characterized at the protein level. It is expressed in very restricted regions of the adult brain. From birth onwards, S100A5 remained a neuronal-specific protein, only located in a subpopulation of neurons in the spiral ganglion.

配方:PBS

分子量:12.3 kDa

序列:Met 1-Lys 93

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.