
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Transferrin Receptor / TFRC / CD71 Protein (His tag)

产品货号: TP07639

第三版

别名:CD71;p90;T9;TFR;TFR1;TR;TRFR

描述:Transferrin receptor protein 1, also known as transferrin receptor, Trfr, p90, CD71 and TFRC, is a single-pass type II membrane protein which belongs to the peptidase M28 family and M28B subfamily. TFRC / CD71 is a membrane-bound protein expressed in larger amounts in proliferating. The specific expression of TFRC can represent a diagnostic tool or a therapeutic target in solid tumours expressing this antigen. Transferrin receptor is necessary for development of erythrocytes and the nervous system. TFRC / CD71 is regulated by cellular iron levels through binding of the iron regulatory proteins, IRP1 and IRP2, to iron-responsive elements in the 3'-UTR. Up-regulated upon mitogenic stimulation. TFRC / CD71 represents a marker of malignant transformation in the pancreas that could be applied as potential diagnostic and therapeutic target.

配方:PBS

分子量:77.4 kDa

序列:Cys 89-Phe 760

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.