
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Glypican 5 / GPC5 Protein (His tag)

产品货号: TP07645

第三版

别名:GPC5;bA93M14.1;glypican proteoglycan 5;glypican-5

描述:Glypican-5 (GPC5), is a cell membrane protein which belongs to the glypican family. The glypicans compose a family of glycosylphosphatidylinositol-anchored heparan sulfate proteoglycans that may play a role in the control of cell division and growth regulation. So far, six members (Glypican-1/GPC1, Glypican-2/GPC2, Glypican-3/GPC3, Glypican-4/GPC4, Glypican-5/GPC5, Glypican-6/GPC6) of this family are known in vertebrates. In adult, Glypican-5 is primarily expressed in the brain. It is also detected in fetal brain, lung and liver. Glypican-5 enhances the intracellular signaling of FGF2 and HGF. It alters the cellular distribution of FGF2. The properties of Glypican-5 make it an attractive target for therapeutic intervention in rhabdomyosarcomas and other tumors that amplify and/or overexpress its gene. Glypican-5 is over-expressed in lymphoma cell lines that had shown amplification. It is a likely target for amplification, and that over-expression of GPC5 may contribute to development and/or progression of lymphomas and other tumors.

配方:PBS

分子量:60.5 kDa

序列:Met 1-Thr 554

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.