



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse REG2 / REG-2 Protein (Fc tag)

产品货号: TP07685

第三版

别名:REG2;RegII

描述:The REG2 Anticoagulation System consists of pegnivacogin, a subcutaneously administered aptamer factor IXa inhibitor, and its intravenous active control agent, anivamersen. The first-in-human study of the effect of the REG2 Anticoagulation System on thrombin generation demonstrates concentration-dependent suppression of thrombin generation that is reversible without rebound effect, as measured by the CAT assay. Full Name

配方:PBS

分子量:43.6 kDa

序列:Met1-Ala173

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.