

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse CRELD1 Protein (His tag)

产品货号: TP07691

第三版

别名: AI843811

描述: CRELD1 is a transmembrane glycoprotein. Epidermal growth factor (EGF) like domain exists in CRELD1. EGF-like repeats are a class of cysteine-rich domains that mediate interactions between proteins of diverse function. EGF domains are found in proteins that are either completely secreted or have transmembrane regions that tether the protein to the cell surface. CRELD1 contains a 333 amino acid (aa) extracellular domain (ECD), two tandem transmembrane segments, and a second ECD of 15 aa. Defects in CRELD1 may cause susceptibility to atrioventricular septal defect type 2 which results in a persistent common atrioventricular canal.

配方: PBS

分子量: 37.7 kDa

序列: Met1-Glu362

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.