

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human CGA / HCG Protein

产品货号: TP07716

第三版

别名:CG-ALPHA;FSH alpha;FSHA;GPHa;GPHA1;HCG;LHA;TSHA

描述:The four human glycoprotein hormones chorionic gonadotropin (CG), luteinizing hormone (LH), follicle stimulating hormone (FSH), and thyroid stimulating hormone (TSH) are dimers consisting of alpha and beta subunits that are associated noncovalently. The alpha subunits of these hormones are identical, however, their beta chains are unique and confer biological specificity. The protein encoded by this gene is the alpha subunit and belongs to the glycoprotein hormones alpha chain family. Two transcript variants encoding different isoforms have been found for this gene.

配方:PBS

分子量:10.2 kDa

序列:Met 1-Ser116

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.