

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Recombinant Human Cystatin A / CSTA Protein (His tag)

产品货号: TP07747

第三版

别名:AREI;STF1;STFA

**描述:**Cystatin-A, also known as Cystatin-AS, Stefin-A and CSTA, is a cytoplasm protein which belongs to the cystatin family. Cystatin-A / CSTA is a cysteine proteinase inhibitor with a molecular mass of 11 kDa, and is located mainly in the keratohyaline granules of the stratum granulosum and the cornified envelope of the stratum corneum in the epidermis. The cystatins are a family of cysteine protease inhibitors with homology to chicken cystatin. Cystatins are physiological inhibitors of cysteine proteinases which are widely distributed in human tissues and fluids. Cystatins typically comprise about 115 amino acids, are largely acidic, contain four conserved cysteine residues known to form two disulfide bonds. Cystatins may be glycosylated and / or phosphorylated, with similarity to fetuins, kininogens, stefins, histidine-rich glycoproteins and cystatin-related proteins. Some of the members are active cysteine protease inhibitors, while others have lost or perhaps never acquired inhibitory activity. Cystatins mainly inhibit peptidases belonging to peptidase families C1 (papain family) and C13 (legumain family).

配方:PBS

分子量:12.4 kDa

序列:Ile 2-Phe 98

纯度:&gt; 95% by HPLC

浓度:

内毒素:&lt;1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.