

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Frizzled-4 / FZD4 / CD344 Protein (Fc tag)

产品货号: TP07791

第三版

别名:CD344;EVR1;FEVR;Fz-4;Fz4;FZD4S;FzE4;GPCR;hFz4

描述:This gene is a member of the frizzled gene family. Members of this family encode seven-transmembrane domain proteins that are receptors for the Wingless type MMTV integration site family of signaling proteins. Most frizzled receptors are coupled to the beta-catenin canonical signaling pathway. This protein may play a role as a positive regulator of the Wingless type MMTV integration site signaling pathway. A transcript variant retaining intronic sequence and encoding a shorter isoform has been described, however, its expression is not supported by other experimental evidence.

配方:PBS

分子量:43.3 kDa

序列:Met 1-Glu180

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.