
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human COL9A1 Protein (Fc tag)

产品货号: TP07817

第三版

别名:DJ149L1.1.2;EDM6;MED;STL4

描述: Ubiquinone biosynthesis protein COQ7 homolog, also known as Coenzyme Q biosynthesis protein 7 homolog, Timing protein clk-1 homolog and COQ7, is a mitochondrion inner membrane and peripheral membrane protein which belongs to the COQ7 family. It is expressed dominantly in heart and skeletal muscle. COQ7 is synthesized as a preprotein that is imported into the mitochondrial matrix, where the sequence is cleaved off and the mature protein becomes loosely associated with the inner membrane. This enzyme is responsible for the hydroxylation of 5-demethoxyubiquinone to 5-hydroxyubiquinone. Human COQ7 protein is mostly helical, and contains an alpha-helical membrane insertion. It has a potential N-glycosylation site, a phosphorylation site for protein kinase C and another for casein kinase II, and three N-myristoylation sites. COQ7 is involved in lifespan determination in ubiquinone-independent manner. It is also involved in ubiquinone biosynthesis. COQ7 is potential central metabolic regulator.

配方: PBS**分子量:** 59.2 kDa**序列:** Met 1-Pro328**纯度:** > 95% by HPLC**浓度:****内毒素:** <1.0 EU per 1 ug of protein (determined by LAL method)**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.