

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Recombinant Human A1BG / alpha 1B-Glycoprotein Protein (His Tag)

产品货号: TP07844

第三版

别名:A1B;ABG;GAB;HYST2477

**描述:**Alpha-1B-glycoprotein (A1BG) is an approximately 80 kDa molecule that circulates in the serum. Mature human A1BG consists of five tandem Ig like domains and shares 44% amino acid sequence identity with mouse and rat A1BG. Alternative splicing generates an isoform that lacks the signal peptide and the first Ig like domain. A1BG binds with high affinity to the cysteine rich protein CRISP 3. In pancreatic carcinoma, it is elevated in pancreatic juice as well as in the serum relative to normal patients. A1BG is abnormally excreted into the urine during idiopathic nephrotic syndrome (INS). Multiple breakdown fragments of A1BG are also present in INS urine, including one of 14 kDa that is found in steroid resistant but not steroid sensitive cases. A1BG is also elevated in the urine of patients with interstitial cystitis/painful bladder syndrome (IC/PBS).

**配方:**PBS

**分子量:**53.4 kDa

**序列:**Met 1-Ser495

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.