

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human GYPA / CD235a / Glycophorin A Protein (His tag)

产品货号: TP07849

第三版

别名:CD235a;GPA;GPErik;GPSAT;HGpMiV;HGpMiXI;HGpSta(C);MN;MNS;PAS-2

**描述:**Glycophorins A (GYPA) and B (GYPB) are major sialoglycoproteins of the human erythrocyte membrane which bear the antigenic determinants for the MN and Ss blood groups. In addition to the M or N and S or s antigens that commonly occur in all populations, about 40 related variant phenotypes have been identified. These variants include all the variants of the Miltenberger complex and several isoforms of Sta, as well as Dantu, Sat, He, Mg, and deletion variants Ena, S-s-U- and Mk. Most of the variants are the result of gene recombinations between GYPA and GYPB.

**配方:**PBS

**分子量:**9.4 kDa

**序列:**Met 1-Glu91

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.