

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human LILRB1 Protein (His Tag)

产品货号: TP07852

第三版

别名:CD85j;ILT-2;ILT2;LIR-1;LIR1;MIR-7;MIR7

**描述:**Leukocyte immunoglobulin-like receptors (LILRs), also called CD85s, ILTs, or LIRs, are important mediators of immune activation and tolerance that contain tandem immunoglobulin (Ig)-like folds. LILRB1 is encoded within the leukocyte receptor complex on 19q13.4, previously implicated to be a susceptibility region to systemic lupus erythematosus (SLE).LILRB1 binds to a wide range of major histocompatibility complex class I molecules (MHCIs) and transduces negative signals that can, for example, prevent killing of MHC1-expressing cells.LILRB1 transcription is correlated with surface protein expression on NK cells.

**配方:**PBS

**分子量:**48.6 kDa

**序列:**Met 1-His458

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.