
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human CLEC1A / CLEC-1 Protein (Fc Tag)

产品货号: TP07881

第三版

别名:CLEC-1;CLEC1

描述:CLEC1A, also known as CLEC-1, is a member of the C-type lectin/C-type lectin-like domain (CTL/CTLD) superfamily. Members of this family share a common protein fold and have diverse functions, such as cell adhesion, cell-cell signalling, glycoprotein turnover, and roles in inflammation and immune response. CLEC1A contains 1 C-type lectin domain and is expressed preferentially in dendritic cells. It may play a role in regulating dendritic cell function. CLEC1A gene is closely linked to other CTL/CTLD superfamily members on chromosome 12p13 in the natural killer gene complex region. Alternative splice variants have been described but their biological nature has not been determined. CLEC1A is found to be not only expressed in dendritic cells, but also in endothelial cells and in the latter aspect resembles the LOX-1 gene.

配方:PBS

分子量:51.9 kDa

序列:Gln22-Asp280

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.