

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human LSAMP Protein (His Tag)

产品货号: TP07900

第三版

别名:FLJ34254;FLJ35396;FLJ37216;FLJ54658;IGLON3;LAMP

描述:The limbic system-associated membrane protein (LAMP) is a cell surface glycoprotein expressed by cortical and subcortical regions of the mammalian CNS that comprise or receive direct projections from limbic system structures. The 64-68-kDa glycoprotein limbic system-associated membrane protein (LsAMP) is expressed on the surface of somata and proximal dendrites of neurons. These areas perform cognitive and autonomic functions, also learning and memory. The functional analysis indicates that LsAMP acts as a selective adhesion molecule, serving as a guidance cue for specific patterns of connectivity, which underlies the normal development of the limbic system. In animal studies there have been found that rats with increased level of anxiety had 1.6-fold higher expression of LsAMP gene in the periaqueductal gray compared to rats with low level of anxiety, indicating a possible role of LsAMP in the regulation of anxiety.

配方:PBS

分子量:33.4 kDa

序列:Met 1-Asn 315

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.