

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Recombinant Mouse TLR3 / CD283 Protein (His tag)

产品货号: TP07908

第三版

别名: AI957183; Tlr3

**描述:** Toll-like receptor 3 (TLR3) also known as CD283 (cluster of differentiation 283) is a member of the Toll-like receptor family of pattern recognition receptors of the innate immune system. TLR3/CD283 plays a fundamental role in pathogen recognition and activation of innate immunity. TLR3 is a nucleotide-sensing TLR which is activated by double-stranded RNA, a sign of viral infection. TLRs are highly conserved from Drosophila to humans and share structural and functional similarities. They recognize pathogen-associated molecular patterns (PAMPs) that are expressed on infectious agents, and mediate the production of cytokines necessary for the development of effective immunity. The various TLRs exhibit different patterns of expression. This receptor is most abundantly expressed in placenta and pancreas, and is restricted to the dendritic subpopulation of the leukocytes. It recognizes dsRNA associated with viral infection, and induces the activation of NF-kappaB and the production of type I interferons. It may thus play a role in host defense against viruses.

**配方:** PBS**分子量:** 78.7 kDa**序列:** Met 1-Leu 705**纯度:** > 95% by HPLC**浓度:****内毒素:** <1.0 EU per 1 ug of protein (determined by LAL method)**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.