
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Pancreasin / Marapsin / PRSS27 Protein (His tag)

产品货号: TP07915

第三版

别名: CAPH2; MPN; UNQ1884/PRO4327

描述: The name "Pancreasin" because it is transcribed strongly in the pancreas. This secreted, tryptic serine protease, also known as Marapsin or PRSS27 (Protease, serine, 27), which is a member of the peptidase S1 family. Pancreasin is inhibited by benzamidine and leupeptin but resists several classic inhibitors of trypsin. Marapsin was constitutively expressed in nonkeratinizing stratified squamous epithelia of human esophagus, tonsil, cervix, larynx, and cornea. In fact, marapsin was the second most strongly up-regulated protease in psoriatic lesions, where expression was localized to the upper region of the hyperplastic epidermis. Similarly, in the hyperproliferative epithelium of regenerating murine skin wounds, marapsin localized to the suprabasal layers, where keratinocytes undergo squamous differentiation. Marapsin's restricted expression, localization, and cytokine-inducible expression suggest a role in the terminal differentiation of keratinocytes in hyperproliferating squamous epithelia.

配方: PBS

分子量: 30.9 kDa

序列: Met 1-Lys 290

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.